

## 10. Department of Biochemistry and Cell biology

- 1) Zhao, Y., Ishigami, M., Nagao, K., Hanada, K., Kono, N., Arai, H., Matsuo, M., Kioka, N., Ueda, K.: ABCB4 exports phosphatidylcholine in a sphingomyelin-dependent manner. *J. Lipid Res.*, 56, 644-652, 2015
- 2) Kang, H.-N., Xu, M., Rodríguez, V.P., Mefed, K., Hanada, K., Ahn, K.-S., Gangakhedkar, S.J., Pakzad, S.R., Prawahju, E.I., Lee, N., Phumiamorn, S., Nemeč, M., Meng, S., Knezević, I.: Review of the current use and evaluation of cell substrates for producing biologicals in selected countries. *Biologicals*, 43, 153-157, 2015
- 3) Sawada, K., Tomoya T., Hamajima, H., Jayakody, L.N., Hirata, M., Yamashiro, M., Tajima, M., Mitsutake, S., Nagao, K., Tsuge, K., Abe, F., Hanada, K., and Kitagaki, H.: Glucosylceramide contained in mold-cultured cereal confers membrane and flavor modification and stress tolerance to yeast during co-culture fermentation, *Appl. Environ. Microbiol.*, 81, 3688-3698, 2015
- 4) Tanida, I., Shirasago, Y., Suzuki, R., Abe, R., Wakita, T., Hanada, K., Fukasawa, M.: Inhibitory Effects of Caffeic Acid, a Coffee-Related Organic Acid, on the Propagation of Hepatitis C Virus. *Jpn. J. Infect. Dis.*, 68, 268-275, 2015
- 5) Fukasawa, M., Nagase, S., Shirasago, Y., Iida, M., Yamashita, M., Endo, K., Yagi, K., Suzuki, T., Wakita, T., Hanada, K., Kuniyasu, H., Kondoh, M.: Monoclonal antibodies against extracellular domains of claudin-1 block Hepatitis C virus infection in a mouse model. *J. Virol.*, 89, 4866-4879, 2015
- 6) Yamashita, M., Iida, M., Tada, M., Shirasago, Y., Fukasawa, M., Nagase, S., Watari, A., Ishii-Watabe, A., Yagi, K., Kondoh, M.: Discovery of anti-claudin-1 antibodies as candidate therapeutics against hepatitis C virus. *J. Pharmacol. Exp. Ther.*, 353, 112-118, 2015
- 7) Nakajima, M., Nagase, S., Iida, M., Takeda, S., Yamashita, M., Watari, A., Shirasago, Y., Fukasawa, M., Takeda, H., Sawasaki, T., Yagi K., Kondoh, M.: Claudin-1 binder enhances epidermal permeability in a human keratinocyte model. *J. Pharmacol. Exp. Ther.*, 354, 440-447, 2015
- 8) Shirasago, Y., Shimizu, Y., Tanida, I., Suzuki, T., Suzuki, R., Sugiyama, K., Wakita, T., Hanada, K., Yagi, K., Kondoh, M., Fukasawa, M.: Occludin-Knockout Human Hepatic Huh7.5.1-8-Derived Cells Are Completely Resistant to Hepatitis C Virus Infection. *Biol. Pharm. Bull.*, in press
- 9) Hanada, K., Sugiki, T.: *In vitro* assay to extract specific lipid types from phospholipid membranes using lipid-transfer proteins: a lesson from the ceramide transport protein CERT, *Neuromethods*, in press